COLLEGE PARK
GUARDIAN™

friction brake knee

joelho com freio de atrito · Reibbremse Knie · rodilla de freno de fricción
ركبة الكبح الاحتكاكي · fonction frein du genou · travão de atrito para joelho
колесо с фрикционным тормозом · 摩擦制动膝关节

college park


TECHNOLOGY for the HUMAN RACE

technical instructions

instruções técnicas · technische anleitung
instrucciones técnicas · التعليمات الفنية
Instructions techniques · Instruções técnicas
Технические инструкции · 技术说明书

TABLE OF CONTENTS

EN	Technical Instructions	1 - 7
BR	Instruções Técnicas	8 - 14
DE	Technische Anleitung	15 - 21
ES	Instrucciones Técnicas	22 - 28
EU	التعليمات الفنية	29 - 35
FR	Instructions Techniques	36 - 42
PT	Instruções Técnicas	43 - 49
RU	Технические инструкции	50 - 56
ZH	技术说明书.....	57 - 63

KEY COMPONENTS

1. Pyramid
2. Upper Assembly
3. Lower Assembly
4. 30mm Pylon Receiver
5. 4mm Clamp Screw (Torque 10 N-m)
6. Release Lever
7. Locking Arm
8. Release Screw
9. External Lock Kit

PACKAGE CONTENTS


- (1) College Park Guardian Knee
- (1) Knee External Lock Kit

TOOLS RECOMMENDED

- (1) 4mm Hex Key

CLEARANCE

1	Overall height	5.0 in (12.8 cm)
2	Dome to knee center	0.9 in (2.3 cm)
3	Dome to tube end contact	2.5 in (6.3 cm)

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
ADJUSTMENTS

All adjustments can be made using a 4mm Allen wrench.

Note: Patient must be sitting when adjustments are made.

STANCE PHASE

Load Dependence

Controls the brake sensitivity. When load is applied, the knee will not bend until the load is displaced. Turning the load adjustment clockwise will increase the amount of load required to initiate braking.

Factory setting = Screw is at minimum (most brake sensitivity).

Note: The factory setting for load dependence should generally not require adjustment. However, it may be necessary in cases such as patients with higher body weights.

Stance Adjustment

Controls the amount of stance flexion. Adjustments can be made to alter the amount of stance flexion before the brake is activated.

Factory setting = Turn screw clockwise to the point where resistance is felt, then back out about 1/8 turn (counterclockwise).

Note: Stance adjustment works together with load dependence to control the brake function. If decreasing the stance adjustment, it may be necessary to also increase the load dependence adjustment.

⚠ Caution: Do not over tighten. May lead to brake sticking which can cause unstable gait.

SWING PHASE

Extension Assist

Controls the rate of extension assist.

Factory setting = Screw is at midpoint of range.

Note: Knee must be fully flexed to access extension adjustment.



DYNAMIC ADJUSTMENTS

Note: The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

⚠ Caution: Flexion and extension must be possible at all settings.

STANCE PHASE

Load Dependence



SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Brake function too strong; Swing too difficult to initiate	Increase Load Dependence (brake less sensitive)	Turn L clockwise	
Brake function not sufficient; Knee too unstable during stance	Decrease Load Dependence (brake more sensitive)	Turn L counterclockwise	

Stance Adjustment

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Not enough flexion	Increase Stance Flexion	Turn S counterclockwise	
Too much flexion	Decrease Stance Flexion	Turn S clockwise	

SWING PHASE

Extension Assist

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Extension too slow or Excessive heel rise	Increase Extension Assist	Turn E clockwise	
Extension too fast or Heel rise not sufficient	Decrease Extension Assist	Turn E counterclockwise	

EXTERNAL LOCKING OPTION

The College Park Guardian can be converted to a locking knee joint, allowing the patient to control when the knee is locked or unlocked. To achieve this, the lock must be enabled and the external lock kit must be installed. Refer to External Lock Kit Fabrication Instructions for detailed information and mounting instructions.


Stance Control Mode	Lock is disabled; Knee utilizes the brake mechanism during stance
Locking Mode	Lock is enabled; Allows for manual locking and unlocking of knee

Factory Setting = Stance Control Mode

ENABLING LOCKING MODE

Note: Do not remove the release screw or apply Loctite

1. Extend the knee joint.
2. Turn the release screw counterclockwise 1-2 turns. Gently push up on the locking arm until it clicks into place.
3. Apply load to the toe of the prosthesis (slight compression of the knee is required to ensure the lock is engaged).
 - ⚠ **Caution:** Test the knee to ensure it does not flex, prior to patient fitting.
4. Continue with installation of the External Lock Kit.


Raised position

LOCK RELEASE FUNCTION

An extension moment is required to unlock the knee when it's in locking mode. Pulling on the manual level alone will not unlock the knee. The patient must pull upward on the manual lever while simultaneously shifting their weight forward slightly.

DISABLING LOCKING MODE

1. Flex the knee (minimum 5 degrees).
 - ⚠ **Caution:** Forcing the locking arm down without first flexing the knee joint will damage the knee.
1. Turn the release screw clockwise until resistance is felt (do not over-tighten). The locking arm will lower slightly and the knee will now be in stance control mode.
 - ⚠ **Caution:** Test the knee to ensure it swings freely without weight applied, prior to patient fitting.

Lowered position







Visit www.college-park.com/guardian to download remote lock fabrication instructions!

USING THE MANUAL LOCK FEATURE

⚠ Caution: Test the manual lever after installation and all adjustments to ensure proper function.

Note: Instruct the patient on proper use of the manual lever. Make sure they understand the operating positions. Visit www.college-park.com to find fabrication instructions.

		
		
<p>LOCKED POSITION</p> <p>Standing/Walking: The knee will not flex</p>	<p>LOCK RELEASE</p> <p>To Sit Down: Patient should pull upward on the manual lever while shifting their weight forward slightly. This allows the knee to flex. Releasing the lever will return it to the locked position.</p>	<p>TEMPORARY UNLOCK</p> <p>Deactivates the lock function temporarily:</p> <p>To allow Temporary Unlock, remove the tabs from the manual lever slot. The knee will remain unlocked when the lever is in this position and will operate in stance control mode.</p>

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Knee Assembly
- External Lock

Warranty inspection schedule for College Park Guardian: Six months, then annually.

WARNINGS

- Flexion and extension must be possible at all settings.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to: noise, sudden loss of function, brake release sticking, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Do not expose this product to moisture, i.e. fresh water, corrosive materials, salt water or pH extremes. If the knee encounters moisture, wipe it dry using a lint-free cloth.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee brake and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.
- If the locking function is utilized, verify that the manual lever and nylon lanyard are not obstructed by the cosmesis.
- If modifying the prosthesis, ensure that the nylon lanyard is secure to prevent damage from grinding tools.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365


College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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©2019 College Park Industries, Inc. All rights reserved.