

Miami J® Cervical Collar



Features and Benefits

- Sorbatex™ padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps ensure patient compliance
- MRI, CT and X-ray lucent
- Bioengineered to minimise pressure points in key known areas: chin, occiput, trapezius and clavicle
- Comfort adjustment buttons for occiput and sternal pad removable when lying and swallowing

Indications for Use

- Post trauma
- Stable cervical fractures
- Cervical spondylitis
- Herniated cervical disc
- RA and OA cervical spine
- Cervical spondylosis
- Post operative

Ordering Information

Size	Product Code		
	Collar Set*	Collar	Replacement Pads
Super Short**	MJR-200S	MJ-200S	MJP-100
Stout	MJR-200L	MJ-200L	MJP-200L
X-Small	MJR-250	MJ-250	MJP-250
Short	MJR-300	MJ-300	MJP-100
Regular	MJR-400	MJ-400	MJP-100
Tall	MJR-500	MJ-500	MJP-100

*Collar sets come with one additional set of replacement pads.

**Super Short collars are packaged with Regular Back Panels.

Sizing Chart

Speciality sizes: 30% of the population

200S SUPER SHORT/KYPHOTIC 10%	200L & 200B STOUT/BARIATRIC 15%	250 XS 5%	500 TALL 10%
Does your patient have a kyphotic (chin-on-chest) neck?	Does your patient have a very large neck circumference?	Does your patient have a very short, thin neck?	Does your patient have a long, tall neck?
<ul style="list-style-type: none"> • Kyphotic • Chin-on-chest • Geriatric • Osteoporitic and/or ankylosing spondylitic 	<ul style="list-style-type: none"> • Very large neck circumference • Obese or massive shoulders <p><small>*For a neck circumference >25", try the Miami J Bariatric Collar (MJ-200B)</small></p>	<ul style="list-style-type: none"> • Very short, thin neck • Principally female • More prevalent in the Asian population • Rare 	<ul style="list-style-type: none"> • Long, tall neck • Young women • Adolescents

Standard sizes: 60% of the population

400 REGULAR 30%	300 SHORT 30%
Is your patient female?	All other adults
<ul style="list-style-type: none"> • Standard size for mature women • Thin, mature men 	<ul style="list-style-type: none"> • Standard size for mature men • Short-necked women

Phenotype based sizing system makes size selection simple whilst drastically reducing margin of error in fitting a patient. Each collar has bioengineered size specific front and back sections designed to work together for superior immobilisation.

The Miami J front and back panels may be interchanged to accommodate patient anatomy, e.g. a short front may be used with a tall back. Choose collar size based on patient neck height, not patient body size. Super short consist of a regular back panel and a super short front panel.