

Universal PROGLIDE 174

Features and Benefits

The PROGLIDE 174 cervical collar is designed to offer maximum support with the addition of an extended sternal pad for increased occipital support. Featuring the same easy-glide system as the PROGLIDE 172, it aims to be a one-size-fits-all solution.

- Extended sternal pad and occipital support
- Ergonomic adjustable control system
- 12 individual adjustments
- Pulley system cervical collar
- Soft cloth inner liner
- Tracheal opening
- MRI, CT and X-Ray compatible

Indications for Use

- Cervical spine precaution for trauma patients
- Immobilisation for pre and post cervical spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis
- Anterior cervical fusion
- Halo removal
- Dens type 1 cervical fractures of C2 vertebrae
- Anterior discectomy
- Teardrop fracture of the vertebral body
- Cervical strain
- Neural involvement
- Atlanto-occipital dislocation



Ordering Information

Description	Size	Product Code
174 Collar	Universal	PG174
Spare Liner Set	Universal	PGPAD174

Miami J® Select Collar

Features and Benefits

- Easy height adjustability with a tamper-proof locking mechanism - 'Left to Lock'
- Anti-microbial and biocompatible Sorbatex™ padding prevents bacterial build-up
- Universal design
- 7 precise fitting sizes with opportunity to fine-tune
- Miami J® inspired chin tray
- One-piece, anatomically inspired collar back with Flex Edge® designed to reduce occipital pressure
- Easy to don and doff with patented Reproducible Fit option
- Intuitive blue Patient Touch-Points to send patients home with confidence - Blue is You™
- Simplified comfort adjustments with integrated and patented Sternal-Relief Dial
- Large tracheal opening - improved function for ICU patients
- Easy to clean and replace multi-coloured padding
- X-ray and CT lucent, MR-safe
- Compatible with Occian™ Back & Miami JTO® accessories

Indications for Use

- Pre and post c-spine surgery
- C-spine precaution for trauma patients



The Miami J Select is designed to support successful patient outcomes through its wide range of height adjustability.

Featuring a patented, first-to-market locking mechanism to protect against tampering and intuitive blue Patient Touch-Points, this collar encourages patient compliance and confidence.

The anti-microbial Sorbatex padding prevents bacteria build-up while being easy to clean and replace when needed.

Ordering Information

Description	Product Code
Collar Set	MJSR-101
Miami J Select Collar	MJS-101

Accessories	Product Code
Miami J Select Replacement Pads	MJSP-101
Miami J Select Tabs Reproducible Fit Tabs (10 pack)	MJST-101