

X-ROM™ Post-Op Knee Brace



Features and Benefits

- Fast, simple application and adjustment through innovative strap, tab and buckle designs
- Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Four independently telescoping sliders
- Ability to contour to the patient's leg to aid varus and valgus support
- Comfortably accommodates a wide range of patient shapes and sizes (from 4'11" to 6'5")
- Enhanced ankle pads and strategically placed contact points between the leg and brace help reduce slippage
- Anti-microbial pads help neutralise odour while increasing comfort and durability

Indications for Use

- To aid immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs

Ordering Information

Size	Colour	Product Code
Universal	Black	11-2181-9



X-ROM™ Post-Op
 Knee Brace
Page 88