

## Lindhe Xtend Connect Component Build Dimensions

	Combination	Description	Functional Height (A)	Effective Height (B)	Total Height (C)	Width Centre to Outer Wheel Locked (D)	Total Width Wheel Locked (E)
1.	NA100-Q4 & NA100-BF	Xtend Connect Quick Fit unit, 4-hole adapter & Xtend Connect Base Unit, female pyramid adapter	19mm	39mm	39mm	48mm	74mm
2.	NA100-Q4 & NA100-BM	Xtend Connect Quick Fit unit, 4-hole adapter & Xtend Connect Base Unit, male pyramid adapter	19mm	26mm	39mm	48mm	74mm
3.	NA100-Q30 & NA100-BF	Xtend Connect Quick Fit Unit, tube clamp 30mm & Xtend Connect Base Unit, female pyramid adapter	19mm	36mm	66mm	38mm	58mm
4.	NA100-Q30 & NA100-BM	Xtend Connect Quick Fit unit, tube clamp 30mm & Xtend Connect Base Unit, male pyramid adapter	19mm	23mm	66mm	38mm	58mm
5.	NA100-QF & NA100-BF	Xtend Connect Quick Fit unit, female pyramid adapter & Xtend Connect Base Unit, female pyramid adapter	19mm	52mm	52mm	38mm	58mm
6.	NA100-QF & NA100-BM	Xtend Connect Quick Fit unit, female pyramid adapter & Xtend Connect Base Unit, male pyramid adapter	19mm	39mm	52mm	38mm	58mm

