

The Terrain foot was created to achieve anatomical ankle motion and energy return in a low-profile design. College Park use their tri-axial technology expertise in a new way. Versatile and service-free, the Terrain features a carbon fibre base and permanent bumpers inside the ankle housing for a progressive, comfortable response during low to high impact activities.

Waterproof* (fresh water only)

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Tri-axial motion, a unique method achieving full movement on all three planes
- A more neutral feel due to its anatomical ankle position and pivot point
- Supported anatomical ankle motion and flexible platform for uneven ground compliance and a balanced energy return
- Maintenance free
- Waterproof* (fresh water only)



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Ordering Example - NTN-LCEN-27-G6						Weight (kg)	Sizes (cm)	Low/Mod	High
Part ID	Side	Shell Colour	Mounting	Size	Firmness Category				
NTN	L		EN		G	0-65	21-30	1	2
NTN	R		EN		G	65-82		2	3
				21-30cm		82-100		3	4
						101-125		4	5
						126-150		5	6
						151-166	26-30	7	7

Spare Parts Ordering Information

Foot Shell	Size (cm)	Wedge Kit	Size (cm)
NHZ-S ■ ● ▲	21-30	NHZ-HWK-2125	21-25
		NHZ-HWK-2630	26-30
CAP	Size (cm)	CPI Spectra Sock	Size (cm)
NBE-C ■ ● -2526	21-25	NHZ-SS-2125	21-25
NBE-C ■ ● -2730	26-30	NHZ-SS-2630	26-30
		NHZ-SSK-2125	21-25 (6 pack)
		NHZ-SSK-2630	25-26 (6 pack)

Shelltread	Size (cm)
NHZ-T ■ N-▲	21-30
Neutral	

Product Specification

Mounting	Sizes	Weight Limit	Build Height
Endo (IP)	21-30cm	21-25cm	149 kg
		26-30cm	166 kg
Weight of Unit	Impact Level	Warranty	Foot Shell
815 g Size 26 inc. foot shell	M H	3 years foot shell - 6 months	C T B J Caucasian, Tan, Brown or Jet Black

