

Designed for superior comfort and rotation for all types of activities. The Voyager allows users to experience a more natural feel with its vertical shock absorption, built-in torsion adapter and responsive springs. Constructed with carbon fibre and College Park's Intelliweave® composites, the Voyager provides a comfortable and natural walking style.

Waterproof* (fresh water only)

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Integrated torsion adapter gives +/- 20° of rotation
- Vertical shock absorption allows for 10mm of movement, reducing impact and socket pressures
- High level of energy return
- Waterproof* (fresh water only)
- Up to 150kg patient weight limit**



Integrated torsion adapter for +/- 20° rotation

Vertical shock unit



High impact



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Ordering Example - Nvy-LCEN-23-G4

Part ID	Side	Shell Colour	Mounting	Size	Firmness Category
Nvy	L		EN		G
Nvy	R		EN		G

Caucasian C
 Tan T
 Brown B
 Jet Black J

21-30cm

Weight (kg)	Sizes 21-25cm	Sizes 26-30cm
0-65	1	1
65-82	2	2
82-100	3	3
101-125	4	4
126-150	-	5

Spare Parts

Ordering Information

Foot Shell	Size (cm)
NHZ-S ■ ● ▲	21-30

■ = "L" or "R" for Left or Right
 ● = Colour (C = Caucasian, T = Tan, B = Brown)
 ▲ = Foot size in cm (21-30)

Shelltread	Size (cm)
NHZ-T ■ N-▲	21-30
Neutral	

Wedge Kit	Size (cm)
Nvy-HWK-2125	21-25
Nvy-HWK-2630	26-30

CPI Spectra Sock	Size (cm)
NHZ-SS-2125	21-25
NHZ-SS-2630	26-30
NHZ-SSK-2125	21-25 (6 pack)
NHZ-SSK-2630	25-26 (6 pack)

Product Specification

Material	Intelliweave® , Carbon Fibre	Unit Weight	1,302g Size 26 inc. foot shell
Weight Limit**	21 - 25cm: 125kg 26 - 30cm: 150kg	Activity Level	K3-K4
Warranty	3 Years foot shell - 6 months	Impact Level	M H
		Footshell	C T B J Caucasian, Tan, Brown or Jet Black Standard Width

