



An infant's first cosmetic prosthesis can be fitted from as early as two to six months of age.

This lightweight solution of a foam-filled, PVC cosmetic glove can assist child development by aiding sitting, balance, crawling and play activities.

- Upper Limb Prosthetics
- Ordering Information
- Myoelectric Hands
- Myoelectric Elbows
- Myoelectric Power & Inputs
- Cable Operated & Spring Hands
- Functional Devices
- Cosmetic Hands & Gloves
- Modular Arm Systems
- Wrist Units
- Mech. Elbow Units & Joints
- Suspension & Control Systems
- Materials & Tools

- Upper Limb Prosthetics
- Ordering Information
- Myoelectric Hands
- Myoelectric Elbows
- Myoelectric Power & Inputs
- Cable Operated & Spring Hands
- Functional Devices
- Cosmetic Hands & Gloves
- Modular Arm Systems
- Wrist Units
- Mech. Elbow Units & Joints
- Suspension & Control Systems
- Materials & Tools

Infant Foam Filled Gloves Ordering Information

- Available in Standard PVC only
 - 19 colours available - colour swatch E24654 applies
 - The below sizes are only available as Foam Filled Gloves
 - See manufacturing guide for Infant Foam Filled Gloves overleaf
- Note:** # Denotes colour.

Infant Foam Filled Gloves Colour Options

E0	E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12	E13	E14	E15	E16	E17	E18
----	----	----	----	----	----	----	----	----	----	-----	-----	-----	-----	-----	-----	-----	-----	-----

	Part no.	Width Across Knuckles	Palm Circumference	Middle Finger Length	Length of Foam Fill (thumb tip to base of glove)	Left/ Right
1	AH45/60#	38mm (1 1/2")	95mm (1 1/2")	45mm (1 3/4")	152mm (6")	Left
	AH45/61#	38mm (1 1/2")	95mm (3 3/4")	45mm (1 3/4")	152mm (6")	Right
2	AH45/25#	41mm (1 5/8")	99mm (3 7/8")	48mm (1 7/8")	178mm (7")	Left
	AH45/26#	41mm (1 5/8")	99mm (3 7/8")	48mm (1 7/8")	178mm (7")	Right
3	AH45/64#	41mm (1 5/8")	105mm (4 1/8")	35mm (1 3/8")	178mm (7")	Left
	AH45/65#	41mm (1 5/8")	105mm (4 1/8")	35mm (1 3/8")	178mm (7")	Right
4	AH45/13#	45mm (1 3/4")	105mm (4 1/8")	49mm (2")	203mm (8")	Left
	AH45/14#	45mm (1 3/4")	105mm (4 1/8")	49mm (2")	203mm (8")	Right
5	AH45/62#	45mm (1 3/4")	111mm (4 3/8")	45mm (1 3/4")	203mm (8")	Left
	AH45/63#	45mm (1 3/4")	111mm (4 3/8")	45mm (1 3/4")	203mm (8")	Right
6	AH45/27#	48mm (1 7/8")	118mm (4 5/8")	54mm (2 1/8")	203mm (8")	Left
	AH45/28#	48mm (1 7/8")	118mm (4 5/8")	54mm (2 1/8")	203mm (8")	Right
7	AH45/19#	49mm (2")	124mm (4 7/8")	57mm (2 1/4")	210mm (8 1/4")	Left
	AH45/20#	49mm (2")	124mm (4 7/8")	57mm (2 1/4")	210mm (8 1/4")	Right

