



Assistive Technology
KEO

 **steeper**
Creating life's turning points, together

KEO

Lightweight and portable home control system

An easy to use, wireless environmental control system designed to enable independent living. Lightweight and portable, KEO facilitates a huge range of everyday home automation tasks, such as opening curtains, answering telephones and changing TV channels.



Home Control

Facilitating a large range of everyday tasks, the KEO can be configured to operate a number of functions within the user's household, granting independence, control and freedom.



Easy to use

With a range of different configuration

and operational devices to control the KEO with, the user has complete control over their environment. This simple to use device, is basic in design, but allows complex control over the user's home at the touch of a button. With the KEO, there are no complicated instruction manuals, instead the user can take back control over their home quickly and easily.



Key Features

- Operated by direct touch or up to two switch scanning control options.
- Ability to program up to 250 infrared codes to control many different household appliances, such as TVs, DVD players, stereos, loud-speaking telephones, intercoms and computers.
- Compatible with existing Steeper infrared peripherals.
- Use with Bluetooth to allow mobile phone control.