



Assistive Technology
pererro+

 **steeper**
Creating life's turning points, together

pererro+

Easy to use switch and voice control for Apple devices

pererro+ by Steeper is a simple, compact, plug and play adapter that works effortlessly with any Apple appliance. It provides users with access to all features, without the need to touch the screen.



Easy to use

Simply use a switch to scan through the on screen options and select the desired function - it really is that easy.



Communication for everyone

Apple devices are the market leaders in mobile communication and we believe this technology should be available to everyone.

Message, call or FaceTime your friends, share that web page on social media, relax and read a book or combine with environmental control apps like evoassist for ultimate control.



Key features

- Plug and play, works via iOS switch control
- Up to four switch access control of iPad, iPad mini, iPod touch and iPhone
- Lightning connector
- Turn On/Off Voice Over
- No external battery, no Bluetooth battery drain
- Powers directly from your device
- 2 x 3.5mm jack sockets for universal switch input
- 'Charge through' ability to charge your device without disconnecting
- Just 33 x 33 x 9 mm

Tel: +44 (0) 113 207 0449 Email: assist@steepergroup.com www.steepergroup.com

STPAT007